

## Veterinary Size - Empty Porcine Gelatin Capsules Physical Specifications


Subject to revision. Contact us for latest specifications

Capsule Size	Typical Fill Weights (gm) Powder Density			Body Volume (ml)	Locked Length ± 1 (mm)		External Diameter ± 0.1 (mm)	Cut Length ± 0.6 (mm)	Single Wall Thickness ± 0.02 (mm)	Dome Thickness ± 0.05 (mm)	Weight (Avg. of 100) ± 10% (gm)
	0.45 (light)	0.7 (typical)	1.00 (heavy)								
Su07	12.6	19.6	28.0	28.0	87.6	Cap	23.30	30.0	0.22	0.27	2.08
						Body	22.70	76.8	0.22	0.27	
07	10.8	16.8	24.0	24.0	77.2	Cap	23.30	30.0	0.22	0.27	1.92
						Body	22.70	66.5	0.22	0.27	
10	8.1	12.6	18.0	18.0	63.8	Cap	23.30	30.0	0.22	0.27	1.69
						Body	22.70	53.0	0.22	0.27	
11	4.5	7.0	10.0	10.0	48.0	Cap	20.68	26.0	0.22	0.27	1.17
						Body	20.10	38.5	0.22	0.27	
12el	3.4	5.3	7.5	7.5	57.0	Cap	15.27	21.0	0.20	0.24	0.85
						Body	14.67	51.0	0.19	0.24	
12	2.3	3.5	5.0	5.0	40.5	Cap	15.27	21.0	0.20	0.24	0.69
						Body	14.67	33.5	0.22	0.24	
13	1.4	2.2	3.2	3.2	30.0	Cap	15.26	18.0	0.20	0.24	0.51
						Body	14.67	23.0	0.22	0.24	

Actual Fill Weights will vary and depend on powder characteristics



For larger capsules sizes see reverse



Capsule volume in ml reflects the capsule body only and does not include the cap dome.

1 3/4; 1 1/2; 1; 3/8; 1/2; 1/4 and 1/8 oz are common terms for the capsule sizes in USA and are not a representation of capsule volume.







## King Size Capsules


Capsule Size	Typical Fill Weights (gm) Powder Density			Body Volume (ml)	Locked Length ± 1 (mm)		External Diameter ± 0.1 (mm)	Cut Length ± 0.6 (mm)	Single Wall Thickness ± 0.05 (mm)	Dome Thickness ± 0.12 (mm)	Weight (Avg. of 100) ± 12% (gm)
	0.45 (light)	0.7 (typical)	1.00 (heavy)								
32150	50	77	110	110	165	Cap	33.40	38	0.28	0.50	8.34
						Body	32.56	150	0.30	0.55	
35110	41	63	90	90	125	Cap	35.50	42	0.28	0.50	6.50
						Body	34.63	109	0.30	0.55	
23100	16	26	36	36	110	Cap	23.30	30.0	0.24	0.27	2.66
						Body	22.70	100	0.23	0.27	

Actual Fill Weights will vary and depend on powder characteristics

For smaller capsule sizes see reverse



Veterinary Capsule Filling Machine  
13 - Su07



ProFill Vet  
13 - Su07



Vet Clinic Filler  
13 - Su07, 32150, 35110